BRAND

Your Partner in the Lab. Worldwide

► Introduction PLT unit

ISO DIN 8655

EUROPÄISCHE NORM
EUROPEAN STANDARD
NORME EUROPÉENNE

EN ISO 8655-2

September 2002

ICS 17.060

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Ref. Nr.: EN ISO 8655-2:2002 D

How often have measuring instruments to be checked?


- Measuring instruments have to be tested at regular intervals
- We recommend about 3 to 12 month for Liquid Handling instruments.
- According to DIN ISO 10012 (Quality assurance requirements for measuring equipment) the confirmation interval depends on the individual applications.

GLP

GLP guidelines

Maintenance, calibration, testing and validation

Equipment shall be adequately inspected, cleaned, and maintained. Equipment used for generation, measurement, or assessment of data shall be adequately tested, calibrated and/or standardized.

Time interval for calibration, re-validation and testing

The frequency for calibration, re-validation and testing (performance verification) depends on the instrument itself, the recommendations from manufacturers of equipment, laboratory experience, and the extent of use.

but....

Every certificate represents the result of a calibration at a specific time only.

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Did you know?


- Most pipettes, independently from the manufacturer, fail over time
- ~ 95 % of failures are related to leakage of the pipetting system, leading to the aspiration of wrong volumes (to low!)
- Adjusting the pipette does not solve the problem!
- Failures often occur between calibration intervals.

Did you know?

Typical failure rates for air-interface pipettes

Undetected or difficult to detect leakage leads to:

- Questionable results
- Waste of time to identify during gravimetric calibration cause and effect of the failure of the pipette
- Selection of shorter calibration intervals to be on the safe side
- Increasing of costs due to the shorter calibration intervals.

What is the PLT unit?

PLT unit stands for **P**ipette **L**eak **T**esting unit.

The BRAND leak testing instrument is a differential pressure measuring instrument for determining the leakage rate of air-displacement pipettes.

A low pressure (vacuum) is generated using a vacuum pump and the increase of pressure in a certain time interval is measured. The changes of the pressure in the time interval leads to a leakage rate. The leakage rate corresponds to a certain volume.

What is the PLT unit?

The leak tester cannot substitute the gravimetric testing for the monitoring of measuring instruments. The instrument should be used to check pipettes on a daily basis as a safeguard for the periods between calibrations.

With or without?

Test with or without pipette tip?

The test can take place **with or without** the pipette **tip**. We recommend that the test be performed with a mounted, unused tip so that the entire pipette system is tested.

When a leak has been identified, the test can be repeated without a tip to determine whether the leak arises from the tip cone/tip coupling region.

Dynamic or static?


Dynamic or static test?

In the **static test** the pipette button is not pressed during the test procedure, i.e. the piston doesn't move.

With the **dynamic test**, in addition to the detection of static leaks, one can also determine whether a defective piston (contamination, scratches) has caused a leak. During the measurement time, the pipette button must be pressed down completely and evenly 2-3 times.

We recommend the dynamic test.

Operating elements

Adapters

Single channel with tip

Single channel without tip


Multichannel with or without tip

PE filter

Flow chart

Flow chart

Main Menu

Pipette Type

Self Test

Switch OFF

Settings...

select | back

Main Menu

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Self Test

Switch OFF

Settings...

select | back

Self Test

Quick Test

Full Test

select | back


Self Test

Quick Test

Full Test

select | back


Results

- The leakage rate limits are warning limits!
- Warning limits for all pipettes are stored.

•If the warning limit is reached (display shows FAIL, red light flashing) the missing volume is gravimetrically detectable.

•This corresponds to 25 % of tolerance in the ISO 8655-2.


Limit

•The leak volume limit for a pipette can be calculated from the leakage rate Q_L , considering the dead volume and the aspiration characteristics of the pipette.

•When a pipette is clean, mechanically error-free and passes the test with the PLT unit, the instrument is within the ISO tolerances.

Results

- a) The progress bar is completely filled in to the bottom: The pipette is leaktight
- b) The progress bar is filled in from the top only down to the limit value region: We recommend that a gravimetric test be performed.
- c) The progress bar is filled in to a point well above the limit value: The pipette leaks. A gravimetric test should be performed, and if necessary the unit should be sent for repair.

Target customers

“BRAND” Repair Service Companies (~ 100 companies had repair trainings at BRAND)

External Repair Service Companies in general, like TSL, biomedis, TTE etc.
(PL possible)

QC departments calibrating and repairing pipettes

End-users especially in segments like

- Health care

- Pharmaceutical industry

- Food industry

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© 2002 CEN Alle Rechte der Vervielfältigung, gleich in welcher Form und in welchem Verfahren, sind weltweit den nationalen Mitgliedern von CEN vorbehalten.