

Optional Accessories

Printer	
Data printer	YDP20-OCE
Paper rolls, 5 units, each with 40 m	6906937
Ink ribbon cassette	6906918
Data Cable	
RS232 25pin male USB-cable type A, length approx. 1.8 m	YCC01-USBM2
RS232 25pin male 25pin female, length approx. 1.5 m	7357312
RS232 25pin male 9pin female, length approx. 2.0 m	7357314
RS232 25pin male 9pin female, length approx. 0.5 m	6965619
Density Determination	
Density determination kit for balances with a readability of 0.1 mg	YDK03
In-use Dust Cover	
for models with a rectangular weighing pan	6960ED01
for models with a round weighing pan (∅ 150 mm)	6960ED02
Dust Cover	
for models with a readability of 0.1 mg	6960BP08
Weighing Table	
made of wood with stone plate	YWT09
made of stone, with vibration dampening	YWT03
Console	
Wall mounting console	YWT04
General	
Remote display, reflective (for connection to data interface port)	YRD03Z
External rechargeable battery pack	YRB11Z
Ionizing blower for eliminating static electricity	
– 220 – 240 V AC	YIB01-ODR
– 110 – 120 V AC	YIB01-OUR
Stat-Pen anti-static device for eliminating electrostatic charges on samples and containers (100 – 240 V AC, 50 60 Hz)	YSTP01

Calibration Weights	
Calibration for lab balance model 224 – Proof Line knob weight 200 g, OIML class E2, with DAkkS certificate	YCW522-AC-02
Calibration for lab balance model 124 – Proof Line knob weight 100 g, OIML class E2, with DAkkS certificate	YCW512-AC-02
Calibration for lab balance model 64 – Proof Line knob weight 50 g, OIML class E2, with DAkkS certificate	YCW452-AC-02
Calibration for lab balance model 623 – Proof Line knob weight 500 g, OIML class F1, with DAkkS certificate	YCW553-AC-02
Calibration for lab balance model 423; 323 – Proof Line knob weight 200 g, OIML class F1, with DAkkS certificate	YCW523-AC-02
Calibration for lab balance model 153 – Proof Line knob weight 100 g, OIML class F1, with DAkkS certificate	YCW513-AC-02
Calibration for lab balance model 6202 – Proof Line knob weight 5 kg, OIML class F1, with DAkkS certificate	YCW653-AC-02
Calibration for lab balance model 4202; 3202; 2202 – Proof Line knob weight 2 kg, OIML class F1, with DAkkS certificate	YCW623-AC-02
Calibration for lab balance model 822 – Proof Line knob weight 500 g, OIML class F2, with DAkkS certificate	YCW554-AC-02
Calibration for lab balance model 8201; 5201 – Proof Line knob weight 5 kg, OIML class F2, with DAkkS certificate	YCW654-AC-02
Calibration for lab balance model 2201 – Proof Line knob weight 2 kg, OIML class F2, with DAkkS certificate	YCW624-AC-02



YDP20-OCE, Data printer