

# BRAND Inserts 2in1

Genuine innovation for cell culture



High flexibility – to be used as a hanging or standing insert

cellGrade™ plus treated membranes for improved cell attachment and growth

Easy to use and efficient in your cell culture lab

**NEW!**



# BRAND Inserts 2in1 – genuine innovation

## ■ Flexible

Keep using the multiwell plates or dishes you already use in your cell culture lab – and work with the advantage of a combined hanging and standing insert when using the BRAND Insert 2in1.


When used as a hanging insert, the insert fits nearly all multiwell plates. For bigger wells and culture dishes, simply use the same insert as a standing insert.

## ■ Functional

The cellGrade™ plus surface of the membranes ensures excellent cell attachment and growth. The highly transparent membranes give you good visibility of your cell cultures, and the high pore density contributes to good permeability.


## ■ Easy to use and efficient

Innovative packaging provides for efficient unpacking and reduces waste.



## Easy-to-use packaging:

① Simply open the entire blister by pulling the sealing paper in one diagonal move (lilac arrow).



② A perforation in the sealing paper also allows to open only one line of inserts and leave the others securely sealed (blue arrow).

BRAND Insert 2in1 is available for 6-, 12- and 24-well sizes and in two pore sizes (0.4 µm and 8.0 µm) with PC-membrane or PET-membrane.

Individually packed in single blisters or in a plate-style blister.



BRAND® and the BRAND word and figurative mark are registered trademarks of BRAND GMBH + CO KG, Germany.

Our technical literature is intended to inform and advise our customers. However, the validity of general empirical values, and of results obtained under test conditions, for specific applications depends on many factors beyond our control. Please appreciate, therefore, that no claims can be derived from our advice. The user is responsible for checking the appropriateness of the product for any particular application.

California Residents: For more information concerning California Proposition 65, please go to [www.brand.de/calprop65](http://www.brand.de/calprop65)

Subject to technical modification without notice. Errors excepted.

BRAND GMBH + CO KG · P.O. Box 1155 · 97861 Wertheim · Germany

Tel.: +49 9342 808-0 · Fax: +49 9342 808-98000 · E-Mail: [info@brand.de](mailto:info@brand.de) · Internet: [www.brand.de](http://www.brand.de)

