

# Compact Shaker KS 15 A



## Description

The reliable Orbital Shaker KS 15 A for general shaking tasks for daily use. The loading capacity is 15 kg on a shaking platform of 400 x 300 mm.

### Basic Equipment

Basic device without shaking plate, **without rack system**

### Loading Capacity

Erlenmeyer flasks 100 ml	20 pieces
Erlenmeyer flasks 250 ml	12 pieces
Erlenmeyer flasks 500 ml	6 pieces
Erlenmeyer flasks 1000 ml	4 pieces
Separating funnels 250 ml	1 pieces

You will find further loading possibilities [here](#)

## Details

**Motion:** orbital  
**Shaking platform:** 400 x 300 mm  
**Max. load:** 15 kg  
**Shaking speed:** 30 - 420 rpm, steplessly adjustabl  
**Stroke:** 17 mm  
**Timer:** 0 - 120 min / continuous  
**Order No.:** 6171 000  
**Accessory:** Rack System Combifix KS  
Universal tray KS  
Incubator Hood TH 15

## Technical specifications

**Electrical supply:** 230 V or 115 V, 50 / 60 H please indicate in case of order  
**Enclosure protection:** IP 21  
**Heat emission:** approx. 5 - 10 W  
**Ambient temperature:** 5°C to 50°C  
**Relative humidity:** ~ 85 %  
**Dimensions (w x d x h):** 510 x 490 x 150 mm  
**Weight:** 18 kg