



sartorius

Purchase a new Arium Mini and receive all consumables and installation free of charge

Total value of
promo :
1.170,00 euro



Features & Benefits

Purchase a new arium[®] mini with unique Bagtank technology

For a Daily Consumption < 10 liter

- Simple set-up with plug and play functionality
- Only system with Bagtank technology
- No sanitization required
- Prevents permanent biofilm
- Produce type 1 water from tap or pretreated water

Do you have a
#passionforscience?

Discover it. Share it. Love it.
www.passionforscience.com

arium® mini

arium® mini Ultrapure Water System

Benefit from the most advanced, space-saving system for daily needs of up to 10 liters of ultrapure water.

Avoid time-consuming tank cleaning cycles with hazardous chemicals.

Connect the system to a pretreated feed water source or to a tap to produce Type 1 water.



Water Quality

Type 1 Water

Type 3 Water

arium® mini

arium® mini essential

arium® mini plus

arium® mini
mini essential
mini plus

arium® mini UV/mini essential UV
mini plus UV

Applications

Buffer and media preparation

Histology

Preparation of reagents, blank samples, etc.

ELISA (enzyme-linked immunosorbent assay)

AAS | GF-AAS (atomic absorption spectrometry)

Photometry

SPE (solid phase extraction)

IC (ion chromatography)

LC-MS (liquid chromatography – mass spectrometry)

GS-MS (gas chromatography coupled with mass spectrometry)

HPLC (high-performance liquid chromatography)

Belgium

For further information, please contact:

Voor 't Labo
Kunstdal 21
9900 Eeklo

phone : +32 9 377 36 38
mail : info@voortlabo.be
website : www.voortlabo.be