

2019 promo on single-channel pipettes
Transferpette® S – until December 31, 2019



Special offer for 2019 models

Order your 2019 model Transferpette® S
now at a special price.

Transferpette® S

- Genuine one-hand operation - set and lock volume, pipet, eject
- Distinct precision and effortless operation
- Volume-change protection for safe pipetting

The new Transferpette® S - available at the beginning of 2020:

With a new design, optimized ergonomics and all of the proven excellent features that make Transferpette® S one of the most popular pipettes for demanding applications.

Save now up to
35%*
on each pipette

* Based on recommended sales prices without V.A.T., as per BRAND price list 2019. You will find all details on the reverse side.

BRAND. For lab. For life.



Transferpette® S Promo 2019

Limited until 31/12/2019 · Fax: +49 9342 808-91515 · sales@brand.de

NOW!

Address:

First name	Last name
Institute / Company	
Street	Zip code and City
E-mail	Tel./Fax

Distributor:



Save 35%* on each pipette!

Order Transferpette® S single-channel pipettes during the promotion period and save 35%. Simply enter the desired quantity in the table and send us the completed form by e-mail, fax or mail!

Yes, I would like to participate in the promotion and order Transferpette® S single-channel pipettes at a special price* according to the general terms and conditions of BRAND GMBH + CO KG for orders placed outside shop.brand.de:

Capacity µl	A ≤ ± % µl	CV ≤ % µl	Cat. No.	Quantity	List price €*	Special price €*	Total price €*
Transferpette® S, with DE-M marking, supplied with performance certificate, shelf/rack mount and silicone grease.							
0.1 - 1	2 0.02	1.2 0.012	7047 68		247,50	160,90	
0.1 - 2.5	1.4 0.035	0.7 0.018	7047 69		247,50	160,90	
0.5 - 10	1 0.1	0.5 0.05	7047 70		247,50	160,90	
2 - 20	0.8 0.16	0.4 0.08	7047 72		221,20	143,80	
5 - 50	0.8 0.4	0.3 0.15	7047 73		221,20	143,80	
10 - 100	0.6 0.6	0.2 0.2	7047 74		221,20	143,80	
20 - 200	0.6 1.2	0.2 0.4	7047 78		221,20	143,80	
100 - 1000	0.6 6	0.2 2	7047 80		221,20	143,80	
500 - 5000	0.6 30	0.2 10	7047 82		247,50	160,90	
1000 - 10000	0.6 60	0.2 20	7047 84		247,50	160,90	
Total:							

Date

Signature

Approval

Yes, I agree that BRAND GMBH + CO KG (further referred to as BRAND) may contact me by

e-mail phone mail

for the purpose of product information, news, events, free tickets, satisfaction surveys and to provide product samples (if not applicable, pls. cross out). BRAND is allowed to process my data for the intended purpose. My personal data will not be transferred to third parties.

I can withdraw my consent by short notice at any time without giving any reason and separately for each of the above communication channels. Contact information for my withdrawal is available at the BRAND Privacy Statement on <https://www.brand.de/en/privacy-statement/> or I may withdraw my consent by sending an e-mail to privacy@brand.de.

Date

Signature

* All prices are recommended sales prices without VAT. Savings in percent compared with BRAND price list 2019.

Your order will be delivered through a distributor.

Transferpette®, BRAND®, BRAND. Forlab. For life.® as well as the logo design mark depicted are trademarks of BRAND GMBH + CO KG, Germany. Other brands and trademarks as referenced in the current printed matter are property of their respective third-party owners.

California Residents: For more information concerning California Proposition 65, please refer to www.brand.de/calprop65

Subject to technical modification without notice. Errors excepted.

BRAND GMBH + CO KG · P.O. Box 1155 · 97861 Wertheim · Germany
Tel.: +49 9342 808-0 · Fax: +49 9342 808-98000 · info@brand.de · www.brand.de

